

You have a solution with a concentration $[4.2 \times 10^7 \text{ cells/mL}]$, how would you dilute your cells to get a concentration of $[2.1 \text{ cells/mL}]$.

You dilute 100 μ L of a sample containing $1 \times 10^6 \text{ cells/mL}$ to a total volume of 1.0mL,
a) what is the concentration of the cells in the final (diluted) solution.
b) what is the dilution factor.

You perform a 1/10 dilution of a stock solution (concentration of 10mg/mL), you then take a 1.0mL sample of the dilution, how much solute would there be in the sample solution.

You make a series of five (1/10) dilutions of your stock culture and bacteria, you then plate 1mL onto TSA agar from each of your five dilutions. If you count 100 colonies on your fifth plate, what is the concentration (in cells/mL) of your stock solution.

You have a 20%(v/v) solution of SDS, you want to make 1.0mL of a 1%(v/v) solution of SDS, how much diluent would you add.

If you have a pellet that contains 372.56 nmoles of DNA how would you make a 100 μ M solution.

How would you make 160mL of 1X TBE from a 5X stock.

You need 112ppm calcium in a 3.2L final volume, If you have a stock solution of 3M CaCl₂, how much chlorine do you end up with.

MW (Ca) = 40.08 g/mol

MW (Cl) = 35.43 g/mol

You have a BSA standard that is 3mg/mL, how much of this would you need to add to 6mL of Coomassie Blue to attain a final concentration of [100µg/mL].

18g/L of a particular chemical is required to obtain a 0.2M solution, how much of this chemical is required to obtain a 300µM solution.

How much Tris base is required to make up 600mL of 0.4M Tris. MW (Tris) = 121.1 g/mol.

If 800mg/mL gives you 800ppm how much gives you 400ppm.

If you prepare a 1/5 dilution of a 150mg/mL solution, what is the final concentration of the solution.

You have 10mL of [50mg/mL] Ethidium Bromide, you need to make 5000 μ L of [10mg/mL] working stock, how would you do this.

You need to make 5000uL of [10mg/mL] solution of Ethidium Bromide, how much [50mg/mL] Ethidium Bromide solution would you need.

How would you make 1.2×10^{-3} L of 1X TBE from a 6X stock.

You have a solution with a concentration $[4.2 \times 10^7 \text{ cells/mL}]$, how would you dilute your cells to get a concentration of $[2.1 \text{ cells/mL}]$.

$$C_1 = 4.2 \times 10^7 \text{ cells/mL}$$

$$C_2 = 2.1 \text{ cells/mL}$$

use: $C_1 V_1 = C_2 V_2$

$$\frac{V_1}{V_2} = \frac{C_2}{C_1} = \frac{2.1 \text{ cells/mL}}{4.2 \times 10^7 \text{ cells/mL}}$$

$$\frac{V_1}{V_2} = 5 \times 10^{-8}$$

$$\frac{V_1}{V_2} = \frac{\text{solute volume}}{\text{Final volume}} = \text{dilution}$$

$$DF = \frac{1}{\text{dilution}} = \frac{1}{5 \times 10^{-8}} = 2 \times 10^7$$

use a

$$\text{dilution} = 5 \times 10^{-8}$$

$$DF = 2 \times 10^7$$

You dilute 100 μ L of a sample containing $1 \times 10^6 \text{ cells/mL}$ to a total volume of 1.0 mL,

- a) what is the concentration of the cells in the final (diluted) solution.
- b) what is the dilution factor.

$$C_1 = 1 \times 10^6 \text{ cells/mL}$$

$$V_1 = 100 \mu\text{L}$$

$$V_2 = 1.0 \text{ mL}$$

use: $C_1 V_1 = C_2 V_2$

$$C_2 = \frac{C_1 V_1}{V_2} = \frac{(1 \times 10^6 \text{ cells/mL})(100 \mu\text{L})}{(1.0 \text{ mL})}$$

$$C_2 = 1 \times 10^8 \frac{\text{cells}}{\text{mL}} \frac{\cancel{\mu\text{L}}}{1 \cancel{\text{mL}}} \left(\frac{\cancel{\text{mL}}}{1000 \mu\text{L}} \right)$$

$$C_2 = 1 \times 10^5 \frac{\text{cells}}{\cancel{\mu\text{L}}} \left(\frac{\cancel{\mu\text{L}}}{1000 \mu\text{L}} \right)$$

$$DF = \frac{\text{Final volume}}{\text{solute volume}}$$

$$DF = \frac{1 \cancel{\text{mL}}}{100 \cancel{\mu\text{L}}} \left(\frac{1000 \cancel{\mu\text{L}}}{1 \cancel{\text{mL}}} \right)$$

$$DF = 10$$

$$C_2 = 100 \text{ cells}/\mu\text{L}$$

You perform a 1/10 dilution a stock solution (concentration of 10mg/mL), you then take a 1.0mL sample of the dilution, how much solute would there be in the sample solution.

$$C_1 = 10 \frac{\text{mg}}{\text{mL}}$$

$$\text{dilution} = \frac{1}{10} = 0.1 = \frac{V_1}{V_2}$$

use $C_1 V_1 = C_2 V_2$

$$C_2 = C_1 \frac{V_1}{V_2} = (C_1)(\text{dilution})$$

$$C_2 = (10 \frac{\text{mg}}{\text{mL}})(0.1)$$

$$C_2 = 1 \frac{\text{mg}}{\text{mL}}$$

now ←

use $C = \frac{M}{V}$

$$M = CV = (1 \frac{\text{mg}}{\text{mL}})(1.0 \text{ mL}) = 1 \text{ mg}$$

You make a series of five (1/10) dilutions of your stock culture and bacteria, you then plate 1mL onto TSA agar from each of your five dilutions. If you count 100 colonies on your fifth plate, what is the concentration (in cells/mL) of your stock solution.


$$\text{dilution} = \frac{1}{10} = 0.1$$

$$V_6 = 1 \text{ mL}$$

$$M_6 = 100 \text{ colonies}$$

$$C_6 = \frac{M_6}{V_6} = \frac{100 \text{ colonies}}{1 \text{ mL}}$$

$$C_6 = 100 \frac{\text{colonies}}{\text{mL}}$$



use

$$C_1 V_1 = C_2 V_2$$

$$C_1 = C_2 \left(\frac{V_2}{V_1} \right) = C_2 (\text{dilution}) = C_3 (\text{dilution})^2 = \dots = C_6 (\text{dilution})^5$$

$$C_1 = (100 \frac{\text{colonies}}{\text{mL}}) \left(\frac{1}{0.1} \right)^5 = 1 \times 10^7 \frac{\text{colonies}}{\text{mL}}$$

You have a 20%(v/v) solution of SDS, you want to make 1.0mL of a 1%(v/v) solution of SDS, how much diluent would you add.

$$C_1 = 20\% (v/v) = \frac{20 \text{ mL}}{100 \text{ mL}}$$

use:

$$C_1 V_1 = C_2 V_2$$

$$V_2 = 1.0 \text{ mL}$$

$$C_2 = 1\% (v/v) = \frac{1 \text{ mL}}{100 \text{ mL}}$$

$$V_1 = \frac{C_2 V_2}{C_1}$$

$$V_1 = \frac{\left(1 \frac{\text{mL}}{100 \text{ mL}}\right) (1 \text{ mL})}{\left(20 \frac{\text{mL}}{100 \text{ mL}}\right)}$$

$$\begin{aligned} \text{solvent} &= 1 \text{ mL} - 0.05 \text{ mL} \\ \text{volume} &= 0.95 \text{ mL} \end{aligned}$$

$$V_1 = 0.05 \text{ mL}$$



solute volume

total = 1.0 mL

If you have a pellet that contains 372.56 nmoles of DNA how would you make a 100μM solution.

$$C_1 = 372.56 \frac{\text{n mol}}{\text{pellet}}$$

set: $C_1 = C_2$

$$C_2 = 100 \mu\text{M} = 100 \frac{\mu\text{mol}}{\text{L}}$$

$$\frac{372.56 \text{ n mol/pellet}}{x} = \frac{100 \mu\text{mol}}{1 \text{ L}}$$

$$\begin{aligned} x &= 372.56 \left(\frac{\text{n mol}}{\text{pellet}} \right) \left(\frac{1 \text{ L}}{100 \mu\text{mol}} \right) \\ &= 3.7256 \left(\frac{\text{n mol}}{\text{pellet}} \right) \left(\frac{1 \text{ L}}{10^3 \mu\text{mol}} \right) \\ &= 3.7256 \times 10^{-3} \frac{\text{L}}{\text{pellet}} \left(\frac{10^3 \text{ mL}}{1 \text{ L}} \right) \\ &= 3.7256 \frac{\text{mL}}{\text{pellet}} \end{aligned}$$

$$x = 3.7256 \frac{\text{mL}}{\text{pellet}}$$

How would you make 160mL of 1X TBE from a 5X stock.

$$V_2 = 160 \text{ mL}$$

$$C_2 = 1X$$

$$C_1 = 5X$$

use

$$C_1 V_1 = C_2 V_2$$

$$V_1 = \frac{C_2 V_2}{C_1} = \frac{(1X)(160 \text{ mL})}{(5X)}$$

$$V_1 = 32 \text{ mL}$$

$$\begin{aligned} \text{solvent volume} &= 160 \text{ mL} - 32 \text{ mL} \\ &= 128 \text{ mL} \end{aligned}$$

solute volume

so

You need 112ppm calcium in a 3.2L final volume, If you have a stock solution of 3M CaCl_2 , how much chlorine do you end up with.

MW (Ca) = 40.08 g/mol

MW (Cl) = 35.43 g/mol

$$C_2 = 112 \text{ ppm} = \frac{112 \text{ mg}}{\text{L}}$$

$$V_2 = 3.2 \text{ L}$$

$$C_1 = 3 \text{ M} = \frac{3 \text{ moles}}{\text{L}} \left(40.08 \frac{\text{g}}{\text{mol}} \right) = 120.24 \frac{\text{g}}{\text{L}} \quad \left(\begin{array}{l} \text{this is of Ca} \\ \text{only, because that} \\ \text{is the concentration} \\ \text{of interest.} \end{array} \right)$$

use: $C_1 V_1 = C_2 V_2$

$$\begin{aligned} V_1 &= \frac{C_2 V_2}{C_1} = \frac{\left(\frac{112 \text{ mg}}{\text{L}} \right) (3.2 \text{ L})}{\left(\frac{120.24 \text{ g}}{\text{L}} \right)} = 2.98 \frac{\text{mg}}{\text{g}} \frac{\text{L}}{1-\text{g}} \left(\frac{\text{g}}{1000 \text{ mg}} \right) \\ &= 2.98 \times 10^{-3} \frac{\text{L}}{\text{L}} \left(\frac{1000 \text{ mL}}{\text{L}} \right) = 2.98 \text{ mL (of Ca)} \end{aligned}$$

so

$$\text{CaCl}_2 \Rightarrow V_{\text{Cl}_2} = V_{\text{Ca}} = 2.98 \text{ mL (of Cl}_2\text{)}$$

$$2.98 \text{ mL of Cl}_2$$

You have a BSA standard that is 3mg/mL, how much of this would you need to add to 6mL of Coomassie Blue to attain a final concentration of [100µg/mL].

$$C_1 = 3 \text{ mg/mL}$$

$$V_2 = 6 \text{ mL}$$

$$C_2 = 100 \text{ µg/mL}$$

use:

$$C_1 V_1 = C_2 V_2$$

$$V_1 = \frac{C_2 V_2}{C_1}$$

$$V_1 = \frac{(100 \text{ µg/mL})(6 \text{ mL})}{(3 \text{ mg/mL})}$$

$$V_1 = 200 \frac{\cancel{\text{µg}}}{\cancel{\text{mL}}} \frac{\text{mL}}{1} \frac{\text{mL}}{\text{mg}} \left(\frac{1 \text{ mg}}{1000 \cancel{\text{µg}}} \right)$$

$$V_1 = 0.2 \text{ mL}$$

18g/L of a particular chemical is required to obtain a 0.2M solution, how much of this chemical is required to obtain a 300µM solution.

$$C_{1in} = 18 \text{ g/L}$$

$$C_{1out} = 0.2 \text{ M}$$

$$C_{2out} = 300 \text{ µM}$$

use:

$$\frac{C_{1in}}{C_{1out}} = \frac{C_{2in}}{C_{2out}}$$

$$C_{2in} = C_{2out} \left(\frac{C_{1in}}{C_{1out}} \right)$$

$$C_{2in} = (300 \text{ µM}) \left(\frac{18 \text{ g/L}}{0.2 \text{ M}} \right)$$

$$C_{2in} = 2.7 \times 10^4 \frac{\cancel{\text{µM}}}{1} \frac{\text{g}}{\text{L}} \frac{1}{\cancel{\text{M}}} \left(\frac{\cancel{\text{M}}}{10^6 \cancel{\text{µM}}} \right)$$

$$C_{2in} = 0.027 \frac{\cancel{\text{g}}}{\cancel{\text{L}}} \left(\frac{1000 \text{ mg}}{1 \cancel{\text{g}}} \right)$$

$$C_{2in} = 27 \text{ mg/L}$$

How much Tris base is required to make up 600mL of 0.4M Tris. MW (Tris) = 121.1 g/mol.

$$V_2 = 600 \text{ mL}$$

$$C_2 = 0.4 \text{ M} = \frac{0.4 \text{ mole}}{\text{L}} \left(121.1 \frac{\text{g}}{\text{mole}} \right) = 48.44 \frac{\text{g}}{\text{L}}$$

use

$$C = \frac{M}{V}$$

$$M = CV = \left(48.44 \frac{\text{g}}{\text{L}} \right) \left(600 \text{ mL} \right) \left(\frac{1 \text{ L}}{1000 \text{ mL}} \right)$$

$$M = 29.1 \text{ g}$$

If 800mg/mL gives you 800ppm how much gives you 400ppm.

part A

part B.

Part A:

$$C_1 = 800 \text{ mg/mL}$$

$$C_2 = 800 \text{ ppm}$$

use

$$C_1 V_1 = C_2 V_2$$

$$\text{dilution} = \frac{V_1}{V_2} = \frac{C_2}{C_1}$$

$$= \frac{800 \text{ ppm}}{800 \text{ mg/mL}}$$

$$\text{dilution} = 1 \frac{\text{ppm mL}}{\text{mg}}$$

Part B:

$$C_2 = 400 \text{ ppm}$$

use

$$C_1 V_1 = C_2 V_2$$

$$\text{dilution} = \frac{V_1}{V_2} = \frac{C_2}{C_1}$$

$$C_1 = \frac{C_2}{\text{dilution}} = \frac{400 \text{ ppm}}{1 \frac{\text{ppm mL}}{\text{mg}}}$$

$$C_1 = 400 \frac{\text{ppm}}{1} \frac{\text{mg}}{\text{ppm mL}}$$

$$C_1 = 400 \frac{\text{mg}}{\text{mL}}$$

If you prepare a 1/5 dilution of a 150mg/mL solution, what is the final concentration of the solution.

$$\text{dilution} = \frac{1}{5} = 0.2 = \frac{V_1}{V_2}$$

$$C_1 = 150 \text{ mg/mL}$$

use $C_1 V_1 = C_2 V_2$

$$C_2 = C_1 \left(\frac{V_1}{V_2} \right) = C_1 (\text{dilution})$$

$$C_2 = (150 \text{ mg/mL})(0.2)$$

$$C_2 = 30 \text{ mg/mL}$$

You have 10mL of [50mg/mL] Ethidium Bromide, you need to make 5000 μ L of [10mg/mL] working stock, how would you do this.

$$C_1 = 50 \text{ mg/mL}$$

$$C_2 = 10 \text{ mg/mL}$$

$$V_2 = 5000 \mu\text{L}$$

use =

$$C_1 V_1 = C_2 V_2$$

$$V_1 = \frac{C_2 V_2}{C_1}$$

$$V_1 = \frac{(10 \text{ mg/mL})(5000 \mu\text{L})}{(50 \text{ mg/mL})}$$

$$V_1 = 1000 \frac{\text{mg}}{\text{mL}} \frac{\mu\text{L}}{1} \frac{\text{mL}}{\text{mg}}$$

$$V_1 = 1000 \mu\text{L} \left(\frac{1 \text{ mL}}{1000 \mu\text{L}} \right)$$

take 1mL of the
10mL you have.

$$V_1 = 1 \text{ mL}$$

You need to make 5000 μ L of [10mg/mL] solution of Ethidium Bromide, how much [50mg/mL] Ethidium Bromide solution would you need.

$$V_2 = 5000 \mu\text{L}$$

$$C_2 = 10 \text{ mg/mL}$$

$$C_1 = 50 \text{ mg/mL}$$

use:

$$C_1 V_1 = C_2 V_2$$

$$V_1 = \frac{C_2 V_2}{C_1}$$

$$V_1 = \frac{(10 \text{ mg/mL})(5000 \mu\text{L})}{(50 \text{ mg/mL})}$$

$$V_1 = 1000 \frac{\text{mg}}{\text{mL}} \frac{\mu\text{L}}{1} \frac{\text{mL}}{\text{mg}}$$

$$V_1 = 1000 \mu\text{L} \left(\frac{1 \text{ mL}}{1000 \mu\text{L}} \right)$$

$$\boxed{V_1 = 1 \text{ mL}}$$

How would you make 1.2×10^{-3} L of 1X TBE from a 6X stock.

$$V_2 = 1.2 \times 10^{-3} \text{ L}$$

$$C_2 = 1X$$

$$C_1 = 6X$$

use:

$$C_1 V_1 = C_2 V_2$$

$$V_1 = \frac{C_2 V_2}{C_1}$$

$$V_1 = \frac{(1X)(1.2 \times 10^{-3} \text{ L})}{(6X)}$$

$$V_1 = 2 \times 10^{-4} \text{ L} \left(\frac{1000 \text{ mL}}{1 \text{ L}} \right)$$

$$\boxed{V_1 = 0.2 \text{ mL}}$$

so

solvent volume = $1.2 \times 10^{-3} \text{ L} - 2 \times 10^{-4} \text{ L}$
 $= 1 \times 10^{-3} \text{ L} \left(\frac{1000 \text{ mL}}{1 \text{ L}} \right)$
 $= 1 \text{ mL}$

solute volume